故障APの交換

1.Maintain→Device Replace→「Auto Sync」の状態を確認



旧APのConfigは保存されている事を確認



Network Branch: TS_Japan Site: Lab_Site ~ Impoints Device Replace Impoints Replacement Replacement Records	
Endpoints Device Replace	
Decords	
Sottings O Refresh Auto Sync ① Auto Delete ② TAuto Sync Time	
Maintain Device Name \$ Device Type Device SN \$ Device Model \$ Device Version Device State Last Cfg Upload At Cfg Sync At Actions	Q
Cloud AP 219801A2959199G0000J WA538-JP Release 2449P12 Offline 2023-09-11 11:36:36 2023-09-11 11:36:36 2023-09-11 11:36:36 T	
CLI Helper Total entries: 1 - 1. Page 1 of 1	1
II File System	
II Restore	
II Config Compare	
II Maintenance	
II Health Check	
Device Replace	
Messages	
🖉 System	

3.「Register New Device & Replace」を選択し、「Next」をクリック



4.新APのSNを記入し、「OK」をクリック

	Cloudnet		Network Sn	mart O&M SD-WAN	Service >			Q 🗄 🗸		
0	Network	Branch: TS_Japan Site: Lab_Site 🗸								
Ęò	Endpoints Device Replace		Register New Device	Register New Device						
	Endpoints Monitor Settings Upgrade Private Versions CLI Helper File System Config Compare Maintenance Health Check Tools Device Replace Messages System	Device Replace Replacement Records O Refresh Auto Sync () Auto Delete () (Device Name Device Type WA538-JP Cloud AP Total entries: 1, current entries: 1 - 1. Page 1 of 1	Register New Device Old Device Info: Site: Povice Name: Device Name: Device SN:	Device Model:WA538-JP Device Version:Release 24 Lab_Site WA538-JP 219801A2359199G00010	49P12	× Last 202	t Cfg Upload At 3-09-11 11:36:36	Cfg Sync At 2023-09-11 11:36:36 < 1 >	Actions To Entri Go to	Q

新APに設定を投入する必要は有りません

5.交換を完了



6.「Refresh」をクリックし、最新情報を確認

C	Cloudnet			Network	Smart O&M	SD-WAN Service►			Q (± ,	🗘 🛞 🗐 Manuals	q97196 - - -
۲	Network	Branch: TS_Japan Site: Lab	_Site v								
Ģ	Endpoints	Device Replace									
6	Monitor	Replacement Re	cords								
ŧ¢Ŷ	Settings	🔾 Refresh 🛛 Auto Sync 🗿	🕑 🚺 Auto Delete 🔊	Auto Sync Time	Replace AP						
4	Maintain Upgrade	Device Name 💠	Device Type	Device SN \$	Device Model \$	Device Version	Device State	Last Cfg Upload At	Cfg Sync At	Actions	۹
-		DEVICEREPLACE-OLD-1	Cloud AP	219801A2959199G0000J	WA538-JP	Release 2449P12	Offline	2023-09-11 11:36:36	2023-09-11 11:36:36	₩	
		WA538-JP	Cloud AP	219801A2959199G0001 Q	WA538-JP	Release 2449P12	Online	-		a ∎ ∎	
-		Total entries: 2 , current entries: '	1 - 2. Page 1 of 1						< 1 >	10 Entri v Go to	1
-											
-											
-											
-											
-	Tools										
-	Device Replace										
۰	Messages										
ø	System										

7.「Records」をクリックし、交換結果を確認



8.交換後の動作確認

